

## 164mer Libraries: Testing cfDNA with 164mer Model

### Tanner Baldwin – 08 August 2024

#### Background

Previously, we optimized the R8 recipe to achieve very satisfactory conversion efficiency while causing minimal DNA damage to lambda DNA with pUC19 and 164mer spike-ins. We sent the reagents for this optimized recipe to Dr. Yuval, our Israeli collaborator and they tested the method on cfDNA. Surprisingly, the results showed that our method performed less well than the Zymo DNA Methylation Gold kit when assessed with digital droplet PCR (ddPCR). We do not have a ddPCR machine in Chicago, but we will test our method more rigorously with cfDNA provided by Yushuai as well as a methylated ds164mer spike-in as a control. We will also test our method starting from 5, 0.5, and 0.05 ng cfDNA with dsDNA 164mer as spike in and compare with EM-seq and Zymo gold kit with the same DNA input.

The final round of R8 optimization showed that R8-4 is the optimal recipe in both 55 C for 90 min and 75 C for 20 min. We also found that the reaction with 10% DME achieves best conversion efficiency (although perhaps at the cost of some additional damage). Since R8-4 at 55 C for 90 min showed better performance in terms of both the DNA damage and BS efficiency. We will try this condition. In addition, since our sequencing data in last round showed that adding 10% DME can further improve the BS efficiency, we will also test R8-4 (40 ul) + 5 ul denatured DNA + 5 ul DME to make comparison.

Tanner did the dsDNA ligation for both cfDNA and 164 bp DNA. The QC of the ligated DNA showed the expected pattern, suggesting that the mistake caused by adding ligation did cause serious self-ligation of cfDNA. However, the OCC purification to try to recovery the cfDNA did cause significant cfDNA loss. Based on the QC results, the concentration of the ligated cfDNA is 4 ng/μl and based on calculation the original cfDNA retained is ~3 ng/μl. We decided to do four group of libraries, with each group containing 3 samples starting from 5, 0.5, and 0.05 ng cfDNA.

#### Materials

Name	Ref/Cat Number
NEBNext UltraShear	#E6655
<del>Zymo Methylated and Non-methylated pUC19 DNA Set</del>	<del>D5017</del>
164mer dsDNA oligo (79 ng/μl)	Yiding's Box
KAPA Evo Prep Kit	10096039001
Zymo EZ DNA Methylation Gold Kit (For Desulphonation)	# D5005/6
KAPA HiFi HotStart Uracil+ ReadyMix (For qPCR/PCR)	KK2801/2
AMPure XP Magnetic Beads (For all beads purifications)	Lab Aliquots (Multiple Vendors?)
NEB EM-seq Kit	
Zymo L-Desulphonation Buffer	

## **Design**

### Material Prep Phase

1. Start first with 474 ng of double stranded 164mer oligo
2. Use KAPA EvoPrep Kit to do End Prep and 5mC adaptor ligation, then beads purification eluting with water
3. Separately, start with ~100 ng of Yushuai's cfDNA (4.98 ng/ $\mu$ l)
4. Use KAPA EvoPrep Kit to do End Prep and 5mC adaptor ligation, then beads purification eluting with water
5. Mix the ligated cfDNA and the ligated 164mer spike-in so that the ration is 99.5:0.5. Take out the calculated amount to do denature and then split into 4 groups. The rest is stored at -80 C containing 10 mM Tris buffer 7.5 for future use.

### Treatment Phase

6. Make 10-fold dilutions in series to prepare stock tubes at 5, 0.5, and 0.05 ng/ $\mu$ l concentrations
7. Before doing BS treatment, add 0.5 M KOH to a final conc. of 0.3 M and incubate at 42 C for 20 min to denature all material
8. Split into 5  $\mu$ l samples for treatments adding 45  $\mu$ l of BS recipe
9. Treatments (19 total) – See the sample information table below
10. Using L-desulphonation buffer to conduct desulphonation for all BS treated samples (X total)
11. qPCR and PCR with unique indexes
12. Give to Chang for sequencing

**Sample Information (19 total)**

Sample ID	Recipe ID	Starting Material	Recipe	Time	Temperature
ST	-	0.5 ng	Only denaturing - no BS	-	-
EA	-	5 ng	EM-seq Kit Reagents	~4.5 hours	Varied
EB	-	0.5 ng	EM-seq Kit Reagents	~4.5 hours	Varied
EC	-	0.05 ng	EM-seq Kit Reagents	~4.5 hours	Varied
ZA	-	5 ng	Zymo DNA Methylation Gold Kit	10min/2.5hr	98 C / 64 C
ZB	-	0.5 ng	Zymo DNA Methylation Gold Kit	10min/2.5hr	98 C / 64 C
ZC	-	0.05 ng	Zymo DNA Methylation Gold Kit	10min/2.5hr	98 C / 64 C
5A	R8-4	5 ng	R8-4 (300 µl 68% NH4BS + 6 µl 20M KOH)	90 min	55 C
5B	R8-4	0.5 ng	R8-4 (300 µl 68% NH4BS + 6 µl 20M KOH)	90 min	55 C
5C	R8-4	0.05 ng	R8-4 (300 µl 68% NH4BS + 6 µl 20M KOH)	90 min	55 C
5DA	R8-4	5 ng	R8-4 (300 µl 68% NH4BS + 6 µl 20M KOH)	90 min	55 C
5DB	R8-4	0.5 ng	R8-4 (300 µl 68% NH4BS + 6 µl 20M KOH)	90 min	55 C
5DC	R8-4	0.05 ng	R8-4 (300 µl 68% NH4BS + 6 µl 20M KOH)	90 min	55 C
7A	R8-4	5 ng	R8-4 (300 µl 68% NH4BS + 6 µl 20M KOH)	20 min	75 C
7B	R8-4	0.5 ng	R8-4 (300 µl 68% NH4BS + 6 µl 20M KOH)	20 min	75 C
7C	R8-4	0.05 ng	R8-4 (300 µl 68% NH4BS + 6 µl 20M KOH)	20 min	75 C
7DA	R8-4	5 ng	R8-4 (300 µl 68% NH4BS + 6 µl 20M KOH)	20 min	75 C
7DB	R8-4	0.5 ng	R8-4 (300 µl 68% NH4BS + 6 µl 20M KOH)	20 min	75 C
7DC	R8-4	0.05 ng	R8-4 (300 µl 68% NH4BS + 6 µl 20M KOH)	20 min	75 C

## Ligation of ds164mer (KAPA EvoPrep Kit – Full volume)

---

### 1. End Prep and A-Tailing

To 6 µl of double stranded 164mer DNA (~474 ng) add:

Reagent	Volume per sample (ul)
RNA H <sub>2</sub> O	29 µl
KAPA End Repair and A-Tailing Ready Mix	25 µl

Sample final volume: **60 µl**

Incubate at: **Pre-cool the block to 4C, add sample, then incubate at 35 C for 30 min followed by 65 C for 30 min (Lid @ 80C)**

### 2. Adaptor Ligation

a. To 60 µl of end repaired DNA (~474 ng)

Reagent	Volume per sample (ul)
NEBNext <i>Methylated</i> Adaptor (No Dilution) - <i>NOT the adaptor included in kit</i>	5 µl
KAPA Ligation ReadyMix	10 µl

Sample final volume: **75 µl**

Incubate at: **Pre-cool block to 4 C, add sample, then incubate at 20 C for 5 minutes with lid at 50 C**

b. To **75 µL of sample**, add:

Reagent	Volume per sample (ul)
USER Enzyme	3.0

Sample final volume: **78**

Incubate at: **37 °C for 15 mins with lid at 50 °C**

### 3. AmPure Beads Purification

a. *Warm beads to RT before use! (~30 min)*

i) 0.8x beads ratio – 62.4 µl.

Elute with 94 µl RNA water to collect 93 µl (~5.04 ng/µl)

#### AmPure Beads Purification Procedure

- Vortex beads to resuspend
- Add **X** µl beads to sample (pipette up and down to mix)
- Incubate for up to 5 minutes at RT


- d. Place on magnet and discard supernatant once clear
- e. Add 200  $\mu$ l of fresh 80% ethanol
- f. Incubate at RT for 30 seconds and remove and discard supernatant
- g. Wash again with 200  $\mu$ l 80% EtOH
- h. Air dry the beads for 5 minutes (DO NOT OVER DRY)
- i. Elute with **X**  $\mu$ l RNA H<sub>2</sub>O by incubating at RT for 2 min
- j. Place on magnet and take **X**  $\mu$ l of supernatant to next step

### **Ligation of Yushuai's cfDNA (KAPA EvoPrep Kit – Full volume)**

---

NOTE: Ligation ready mix was mistakenly added in place of ERAT ready mix. OCC was performed by adding 40  $\mu$ l of water to the 60  $\mu$ l of sample, splitting into two, 50  $\mu$ l samples and then following zymo protocol exactly with all brand-new reagents. Elution was done with 18  $\mu$ l and 17.5 was successfully recovered from each sample to give a combined 35  $\mu$ l to repeat the ERAT step using the proper ready mix.

- Zymo OCC claims 16 nt and up will be recovered with 90% efficiency (their data is with ssDNA however), DCC advertises 70-90% recovery of 50ng and up, based on Ampure Beads data recovery would be low and there could be small loss of 160mer even with 2.5x volume



Factoring in adapter dimer, peaks are at 110-130, 300, and 460bp.

Total material with adapter dimer is  $1.1850 \times 4$  (dilution factor) =  $4.74 \times 14.5\text{ul}$  (total volume remaining) = 68.73ng of material total

If we assume the ~130bp adapter dimer comprises ~37% of material (130bp/avg size) then we have 43.299ng of cfDNA after the OCC, ligation, and beads purification.

#### 4. End Prep and A-Tailing

To 35  $\mu\text{l}$  of Yushuia's cfDNA (~160 ng?) add:

Reagent	Volume per sample (ul)
RNA H <sub>2</sub> O	1.5 $\mu\text{l}$
KAPA End Repair and A-Tailing Ready Mix	25 $\mu\text{l}$

Sample final volume: 60  $\mu\text{l}$

Incubate at: **Pre-cool the block to 4C, add sample, then incubate at 35 C for 30 min followed by 65 C for 30 min (Lid @ 80C)**

#### 5. Adaptor Ligation

- To 60  $\mu\text{l}$  of end repaired DNA (~160 ng?)

Reagent	Volume per sample (ul)
NEBNext <i>Methylated</i> Adaptor (No Dilution) - <i>NOT the adaptor included in kit</i>	5 µl
KAPA Ligation ReadyMix	10 µl

Sample final volume: **75 µl**

Incubate at: **Pre-cool block to 4 C, add sample, then incubate at 20 C for 5 minutes with lid at 50 C**

b. To **75 µL of sample**, add:

Reagent	Volume per sample (ul)
USER Enzyme	3.0

Sample final volume: **78**

Incubate at: **37 °C for 15 mins with lid at 50 °C**

## 6. AmPure Beads Purification

a. *Warm beads to RT before use! (~30 min)*

i) 0.8x beads ratio – 62.4 µl.

Elute with 16 µl RNA water to collect 15 µl (~4 ng/µl)

Factoring out the adapter dimer, peaks are at 110-130, 300, and 460bp.

Total material with adapter dimer is  $1.1850 \times 4$  (dilution factor) =  $4.74 \times 14.5\text{ul}$  (total volume remaining) = 68.73ng of material total (including adapter!)

If we assume the ~130bp adapter dimer comprises ~37% of material (130bp/avg size) then we have 43.299ng of cfDNA after the OCC, ligation, and beads purification.

**cfDNA conc. = ~3 ng/µl**

## AmPure Beads Purification Procedure

- k. Vortex beads to resuspend
- l. Add **X** µl beads to sample (pipette up and down to mix)
- m. Incubate for up to 5 minutes at RT
- n. Place on magnet and discard supernatant once clear
- o. Add 200 µl of fresh 80% ethanol
- p. Incubate at RT for 30 seconds and remove and discard supernatant
- q. Wash again with 200 µl 80% EtOH
- r. Air dry the beads for 5 minutes (DO NOT OVER DRY)
- s. Elute with **X** µl RNA H<sub>2</sub>O by incubating at RT for 2 min
- t. Place on magnet and take **X** µl of supernatant to next step

## 7. Split Ligated material into aliquots of 5, 0.5, and 0.05 ng/μl by making a series dilution

- a. Dilute the ligated ds164mer (~5 ng/μl) by taking 1 μl of ds164mer and adding 22 μl water → final conc. ~0.217 ng/μl (not including the adapter)

### Our Samples

<b>Tube 1</b> 5 ng samples desired (~4.69 ng)	14 μl of ligated sample (~42 ng cfDNA) + 3 μl water + 1 μl diluted ds164mer  TOTAL: 18 μl, cfDNA = 42 ng (2.33 ng/μl), ligated ds164mer = 0.217 ng (0.012 ng/μl = ~0.5%)  <b>Take 10 μl for Denaturing, save 8 μl in 10 mM Tris-HCl Buffer pH 7.5 @ -80 C</b> Add 6 μl 20 M KOH to 234 μl RNA water to make 0.5 M KOH Add 15 μl of 0.5 M KOH to the 10 μl of material Incubate in LARGE thermocycler (which has volume >50 μl) at: <b>42 C for 20 minutes with lid at 62 C</b>  TOTAL: 25 μl
<b>Tube 2</b> 0.5 ng samples desired (~0.52 ng)	3 μl Tube 1 + 24 μl water (0.094 ng/μl – 5 μl per sample)  *Take ST sample from this stock*
<b>Tube 3</b> 0.05 ng samples desired (~0.052 ng)	3 μl Tube 2 + 27 μl water (0.0094 ng/μl – 5 μl per sample)

### EM-seq and Zymo Kit Samples

<b>Tube 1</b> 5 ng samples desired (~4.67 ng)	Take 10 μl of the stored material at -80 C (~1.87 ng/μl) <ul style="list-style-type: none"> <li>- Add 10 μl water</li> <li>- Take 12 μl to be the stock for these “5 ng” samples</li> </ul>
<b>Tube 2</b>	Take 2 μl of Tube 1 and add 16 μl water (9-fold dilution) <ul style="list-style-type: none"> <li>-</li> </ul>



0.5 ng samples desired (~ ng)	
<b>Tube 3</b> 0.05 ng samples desired (~ ng)	Take 2 µl of Tube 2 and add 18 µl water (10-fold dilution)

## Homemade BS treatments

### 1. Prepare Reagents

#### 20 M KOH

Reagent	Volume per sample (ul)
KOH	1119 mg
H <sub>2</sub> O	997 ul

Sample final volume: ~997 µl @ ~20 M

#### R8-4 Recipe

Reagent	Volume per sample (ul)
NEW 68% ammonium bisulfite	300
20 M potassium hydroxide	6

Sample final volume: ~306 µl

### 2. Treatments

- For all other treatments, take two tube strips of *easy cut* 200 µl tubes and label them with the sample IDs
- Add 5 µl of denatured material to each sample
  - Note: For DME samples add 5 µl DME here BEFORE adding bisulfite
  - Then in the next step, add only 40 µl of BS recipe (total recipe volume = 45 µl)
- Fill sample tubes with 45 µl of the proper BS recipe
- Treat according to the sample information table
- ~~Preheat the recipes at the proper treatment temperature (as indicated by the sample information table) for up to 2 minutes.~~

- f. ~~Pipette 2.0 µl of starting material (ligated lambda DNA) into each preheated recipe and let the reaction proceed for time specified in the sample information table~~
- g. When the reaction finishes, remove and add the 100 µl of H<sub>2</sub>O that is required for desulphonation. Proceed to desulphonation.

**Notes:**

- DME smells putrid after treating the DNA, worse smell in 75 C than in 55 C
- DME samples are slightly cloudy after adding the BS recipe, after treatment they become clear
- EM-seq DTT reagent had precipitates that cannot be vortexed away

**3. Desulphonation (BS treated samples only) USE DNA KIT – Note any precipitates!**

- a. Prepare X, 1.5 ml tubes and fill them with M-Binding Buffer – 600 µL
- b. Add the 150 µl sample (from treatment step) to the binding buffer and mix well.
- c. Once M-binding buffer and sample have been mixed → Centrifuge Max. speed for 30 s
- d. M-Wash Buffer – 100 µL → Centrifuge Max. speed for 30 s
- e. ~~L-Desulphonation Buffer~~ – 200 µL on column  
M-Desulphonation Buffer was used because L-desulph did not arrive!
- f. Incubate at Rm Temp for 20 mins → Centrifuge Max. speed for 30s
- g. M-Wash Buffer – 200 µL → Centrifuge Max. speed for 30 s
- h. Repeat washing for another 3 times
- i. **Discard → Centrifuge Max. speed for 2 mins**
- j. Elute with H<sub>2</sub>O with 24 µl to get 23 µl for each sample.

**qPCR (19 Samples)**

**1. qPCR**

- a. To **1 µl** of each sample (<1.8 ng), add

Reagent	Volume per 1 sample (ul)	Volume per 20 samples (ul)
KAPA HiFi U+ MM	10	200
Universal primer	1	20

Mixed Index primer	1	20
H <sub>2</sub> O	6	120
<b>20x</b> SYBR Green	1	20
<b>Final volume</b>	<b>19</b>	<b>380</b>

Sample final volume: **20 µl** →

Pre-incubation	95 C; 600 s
X25	95 C; 20 s
	60 C; 20 s
	72 C; 20 s
Melting curves	95 C; 10 s
	65 C; 60 s
	97 C; 1 s

		ST	5A	5B	5C	5DA	5DB	5DC			
			7A	7B	7C	7DA	7DB	7DC			
			ZA	ZB	ZC	EA	EB	EC			

## 2. PCR

a. To **10 µl (< 18 ng)** of sample, add:

Reagent	Volume per sample (µl)	Volume per 24 samples (µl)
---------	------------------------	----------------------------

KAPA HiFi U+ MM	12.5	250
Universal primer	1.25	25
<b>Final Volume</b>	<b>13.75</b>	<b>275</b>
Index primer	1.25	<i>Must be added individually</i>

Sample final volume: **25 ul** →

x cycles # (Calculate based on the qPCR Ct value)	98 C; 45s
	98C; 15 s
	60 C; 30 s
	72 C; 30 s
	72 C; 1 min
	4C

### 3. TWO rounds of AmPure beads purification

- a. 0.6x Beads ratio – 15 ul
- b. H<sub>2</sub>O – 21 ul
- c. 0.6x Beads ratio – 12 ul
- d. H<sub>2</sub>O – 12 ul
  - i. < 2ng/ul: 1ul+2ul H<sub>2</sub>O to QC
  - ii. > 2ng/ul: 0.5+2.5ul H<sub>2</sub>O to QC

9, 13, 15, 17

### Library information


ID	Recipe ID	Recipe	Time	Temperature	qPCR	PCR	Index	Qubit [ng/uL]	QC Conc [ng/uL]
ST	-	Only denaturing - no BS	-	-	8.6	9	B12	3.44	
EA	-	Only denaturing - no BS	~4.5 hours	Varied	9.7	9	C12	2.82	
EB	-	EM-seq Kit Reagents	~4.5 hours	Varied	14.5	16	D12	10.2	
EC	-	EM-seq Kit Reagents	~4.5 hours	Varied	18.4	16*	E12	2.16	
ZA	-	EM-seq Kit Reagents	10min /2.5hr	98 C / 64 C	8.8	9	F12	2.98	
ZB	-	Zymo DNA Methylation Gold Kit	10min /2.5hr	98 C / 64 C	12.5	13	G12	3.8	



<b>ZC</b>	-	Zymo DNA Methylation Gold Kit	10min /2.5hr	98 C / 64 C	16.2	16*	<b>H12</b>	<b>4.54</b>	
<b>5A</b>	R8-4	Zymo DNA Methylation Gold Kit	90 min	55 C	8.4	9	<b>C11</b>	<b>3.5</b>	
<b>5B</b>	R8-4	R8-4 (300 µl 68% NH4BS + 6 µl 20M KOH)	90 min	55 C	12.5	13	<b>D11</b>	<b>3.9</b>	
<b>5C</b>	R8-4	R8-4 (300 µl 68% NH4BS + 6 µl 20M KOH)	90 min	55 C	15.3	16	<b>E11</b>	<b>5.18</b>	
<b>5DA</b>	R8-4	R8-4 (300 µl 68% NH4BS + 6 µl 20M KOH)	90 min	55 C	9.3	9	<b>F11</b>	<b>1.99</b>	
<b>5DB</b>	R8-4	R8-4 (300 µl 68% NH4BS + 6 µl 20M KOH)	90 min	55 C	12.2	13	<b>G11</b>	<b>3.42</b>	
<b>5DC</b>	R8-4	R8-4 (300 µl 68% NH4BS + 6 µl 20M KOH)	90 min	55 C	15.0	16	<b>H11</b>	<b>2.66</b>	
<b>7A</b>	R8-4	R8-4 (300 µl 68% NH4BS + 6 µl 20M KOH)	20 min	75 C	9.0	9	<b>C10</b>	<b>3.7</b>	
<b>7B</b>	R8-4	R8-4 (300 µl 68% NH4BS + 6 µl 20M KOH)	20 min	75 C	12.8	13	<b>D10</b>	<b>3.42</b>	
<b>7C</b>	R8-4	R8-4 (300 µl 68% NH4BS + 6 µl 20M KOH)	20 min	75 C	17.2	16*	<b>E10</b>	<b>3.16</b>	
<b>7DA</b>	R8-4	R8-4 (300 µl 68% NH4BS + 6 µl 20M KOH)	20 min	75 C	10.4	13	<b>F10</b>	<b>6.38</b>	
<b>7DB</b>	R8-4	R8-4 (300 µl 68% NH4BS + 6 µl 20M KOH)	20 min	75 C	13.0	13	<b>G10</b>	<b>1.56</b>	
<b>7DC</b>	R8-4	R8-4 (300 µl 68% NH4BS + 6 µl 20M KOH)	20 min	75 C	16.8	16*	<b>A10</b>	<b>1.91</b>	



All samples that received 9 cycles were incubated too long in the initial denaturing step at 98 C. All samples incubated for

## Post Library QC Results


---











## Discussion

Note: in human cfDNA CHG and CHH motifs are only lightly methylated (~0.5%?)

The post library QC is very encouraging as it reveals the relative peak heights of the three major cfDNA peaks. Our UMBS-seq method shows that longest sized peak (~700) is mostly intact whereas Zymo Gold and EM-seq have depleted this peak size (which is largest and therefore most sensitive to damage).

Sequencing results reveal that the trends are mostly as expected. We see high methylation in promoter regions, and the trend is most prominently revealed by our UMBS-seq method. Also, Em-seq showed a notable drop off in efficiency when moving from 5 ng of input to 0.05 ng. UMBS dropped off only very minimally. This is exactly what we hoped to see as EM-seq is likely less damaging than us, but our chemical method will allow us to stay efficient even at ultra-low inputs.

HOWEVER, something went wrong with the EM-seq kit (maybe the enzymes are expired) because the conversion efficiency was much less than expected. Therefore, we will have to repeat this experiment.